

Micro Manager for Icy



μ-Manager

The Open Source Microscopy Software

Plan

Installation

Configuration

Acquisition

- Snap / Album
- Live
- Multi-D acquisition

Scripting / Development

- Protocols
- (Java) Script
- Plugin

Installation

By default the *μManager for Icy* plugin should be already installed. Verify than you have the last version of the plugin then launch it.

Installation


On first start you should specify the μ Manager's installation folder to Icy.

Be sure you installed a compatible version of μ Manager (currently only version from 1.4.19 to 1.4.23 are supported).



Configuration


As in μ Manager you should then select which configuration to load.



Main window


The main window is very similar to μ Manager with some cleanup and minors changes.

In blue we find the *Configurations Settings* part of μ Manager with *Groups* and *Presets*, in green we have everything about the acquisition itself (actions and settings) and finally in the red part we can find all *μ Manager for Icy* compatibles plugin to extend its features.



Configuration



Configuration can be done exactly as the original μ Manager. It's also possible to access some of the main basic features of μ Manager from the main window menu (by clicking on the top left icon). We can find here for instance loading / saving of settings files, the *Configuration Wizard*, the *Property Browser* and the *Pixel Size Config* of the original μ Manager.



Acquisition

Acquisition part again is very close to the one from μ Manager. We can find here the same actions (in blue), with the same available camera settings (in red) and of course the acquisition information (in green).

The only change are these new parameters (in yellow) which allow to do live and/or basic acquisition directly in 3D.






Snap / Album

As in μ Manager the *Snap* operation create a new image for each acquisition where the *Album* operation will append all acquisitions in the same Sequence. As presented before and unlike the original μ Manager, we can now directly do 3D stack acquisition.


Snap 2D




Snap 3D



Album 2D





Album 3D







Live

Live mode gives you a real time view from the camera as in μ Manager except we can now get a 3D stack view by modifying the parameters.



← LIVE 2D




LIVE 3D →





Live

The advantage of Live 3D is that it can take benefit from the 3D raycasting rendering of VTK to offer a real time 3D view.







Multi-D Acquisition

This plugin correspond to the powerful *Multi-D Acquisition* tool from μ Manager.

The graphical interface is exactly the same as the one we can find in μ Manager except we can now see the acquisition progress.



And other...




Plugin *Remote* is almost the same (except GUI) than the *Stage Position Control* from μ Manager (XY and Z stage position control)



Protocols


Allowing microscope control and image acquisition directly from the protocols !



Protocols - exercise 1

Goal: Design a protocol which can do an image acquisition on 2 channels.


Protocols - exercice 1 - solution



Protocols - exercice 2

Goal: Design a protocol which allow to do a 3D stack acquisition. The stack should contains 10 slices where the first slice is located at Z position = $10\mu\text{m}$ then each slice are separated by $5\mu\text{m}$ space.

Protocols - exercice 2 - solution



Script

µManager for Icy allow to control your microscope from a simple Java script.
WARNING: don't forget to launch the Micro-Manager plugin before !

This simple example will do the following operations :

- Move the XYZ stage at position [5,5,5] µm
- Snap a single image and display it

```
// move the stage to (5, 5, 5)
StageMover.moveXYAbsolute(5, 5)
StageMover.moveZAbsolute(5)

// acquire a single image
image = MicroManager.snapImage()

// create a sequence and display it
sequence = new Sequence(image)
gui.addSequence(sequence)
```

Script / Development

The main classes and methods to know to use the μ Manager API in Icy

MicroManager	classe principale de μ Manager pour Icy
StageMover	classe outil pour gérer le positionnement du microscope

MicroManager.getCore()	Get access to μ Manager core
MicroManager.setExposure(...)	Change camera exposure
MicroManager.snapImage()	Acquire an image and return it
MicroManager.getMetadata()	Retrieve meta data of the last acquired image
MicroManager.startLiveMode()	Start continuous acquisition mode (live)
MicroManager.stopLiveMode()	Stop continuous acquisition mode (live)
MicroManager.startAcquisition(...)	Start multiple acquisition
MicroManager.stopAcquisition(...)	Stop multiple acquisition
MicroManager.getAcquisitionResult()	Retrieve result from multiple acquisition

Script - exercice 1

Goals:

- Move the microscope stage (x,y,z) to [5,5,5]
- Acquire 3 images
- Move the stage of 10 μ m in Z between each acquisition
- Display the result in Icy as a 3D stack image.

Script - exercice 1 - solution

```
z                                =                                0
sequence = new Sequence()                // create the result sequence

StageMover.moveXYAbsolute(5, 5)          // move to position (5, 5, 5)
StageMover.moveZAbsolute(5)

image = MicroManager.snapImage()        // acquire 1 image
sequence.setImage(0, z++, image)        // set it in resulting sequence at position 0
StageMover.moveZRelative(10)            // shift microscope Z position by 10
image = MicroManager.snapImage()        // acquire 1 image
sequence.setImage(0, z++, image)        // set it in resulting sequence at position 1
StageMover.moveZRelative(10)            // shift microscope Z position by 10
image = MicroManager.snapImage()        // acquire 1 image
sequence.setImage(0, z, image)          // set it in resulting sequence at position 2

gui.addSequence(sequence)                // show the sequence in Icy
```

Script - exercice 2

Goals:

- Acquire a Sequence with :
 - 20 frames with exposure time = $10 + (\text{frame_index} * 5)$
 - 10 slices with Z position = $-5 + \text{slice_index}$
- Do an automatic threshold on obtained Sequence and display the result as a ROI (Region Of Interest)

Tips: Automatic threshold can be done using the KMeans method to retrieve the threshold intensity, then we process the threshold itself based on this intensity value (KMeans.computeKMeansThresholds(...) and Thresholder...)

Script - exercice 2 - solution

```
sequence = new Sequence() // create the result sequence
gui.addSequence(sequence) // display it

StageMover.moveXYAbsolute(5, 5) // move microscope to position XY (5, 5)

for(t = 0; t < 20; t++)
{
    MicroManager.setExposure(10 + (t * 5)) // set exposure depending T position

    for(z = 0; z < 10; z++)
    {
        StageMover.moveZAbsolute(-5 + z, true) // set microscope Z position by 10
        image = MicroManager.snapImage() // acquire 1 image
        sequence.setImage(t, z, image) // set it in resulting sequence at position 0
    }
}

value = KMeans.computeKMeansThresholds(sequence, 0, 2, 256) // find threshold value
rois = Thresholder.threshold(sequence, 0, value) // apply threshold and get ROIs

for(i = 0; i < rois.length; i++) // put ROIs on the sequence
    sequence.addROI(rois[i])
```

μManager core access

You can access the internal μManager core and so get access to all the functionalities of the internal μManager API. For instance you can grab value for a specific property and more generally modify some acquisition parameters (see [Programming Guide - Using device properties](#))

μManager core usage (from μManager):

```
core.getProperty(...)
```

μManager core usage (from Icy):

```
MicroManager.getCore().getProperty(...)
```

Script - μ Manager core access

```
importClass(Packages.org.micromanager.utils.MDUtils)

core = MicroManager.getCore()
image = MicroManager.snapImage()
meta = MicroManager.getMetadata()

println("Binning: " + MDUtils.getBinning(meta))
println("Pixel type: " + MDUtils.getPixelType(meta))

bd = core.getProperty("Camera", "BitDepth")
exposure = core.getProperty("Camera", "Exposure")
MicroManager.setExposure(10)
core.setProperty("Camera", "Binning", 2)
```

MicroscopePlugin class

When we develop a new Icy plugin for Micro-Manager it's important to extend the abstract class *MicroscopePlugin* instead of *Plugin* or *PluginActionable*. In this case it's important to respect the following rules:

- Overload the *start()* method (instead of the *run()* method)
- Overload the *shutdown()* method if some specific actions need to be done when plugin is terminated.

Using the *MicroscopePlugin* class assure that μ Manager will be loaded before the plugin starts, also it provides methods as *onSystemConfigurationLoaded()*, *onCorePropertyChanged()* and *onExposureChanged()* to detect configuration changes from μ Manager.

Events

Micro-Manager for Icy adds new events to make life easier for developer.

MicroManager.addAcquisitionListener(...)

Allow to listen acquisition events (start / new image / end).

MicroManager.addLiveListener(...)

Allow to listen events for Live mode (start / new image / end).

StageMover.addListener(...)

Allow to listen events from the microscope stage position (position changed)

So the developer can, for instance, easily start a specific task when receiving a new image during the acquisition.

Plugin - Tutorial project 1

```
public class MyPlugin extends MicroscopePlugin {
    @Override
    public void start()
    {
        try {
            Sequence result = new Sequence();           // Create the resulting sequence
            StageMover.moveXYAbsolute(5, 5);           // Set microscope X and Y positions
            StageMover.moveZAbsolute(5);               // Set microscope Z position
            result.addImage(MicroManager.snapImage()); // Snap an image and add it to result
            StageMover.moveZRelative(10);              // Move the microscope by 10 µm in Z
            result.addImage(MicroManager.snapImage()); // Snap again
            StageMover.moveZRelative(10);              // Move again
            result.addImage(MicroManager.snapImage()); // Then Snap again
            addSequence(result);                       // Finally, show the resulting sequence into Icy
        } catch (Exception e) {
            // Eclipse will ask you to catch the exception, this is caused when we are unable to move the stage
        }
    }
}
```

Plugin - exercice

Objectifs :

- Start the *Live* mode
- Register to receive events from *Live* mode.
- For each received image, display the XY dimension in the output console.

Plugin - exercice - solution

```
public class MyPlugin extends MicroscopePlugin implements LiveListener {
    public void start() {
        try {
            MicroManager.addLiveListener(this); // register listener first
            MicroManager.startLiveMode();      // then start live acquisition
        } catch (Exception e) {
            // we need to catch possible exception here on startLiveMode()
        }
    }

    public void liveImgReceived(IcyBufferedImage image) {
        try {
            JSONObject meta = MicroManager.getMetadata();
            System.out.println("Image size: " + MDUtils.getHeight(meta()) + " x " + MDUtils.getWidth(meta));
        } catch (JSONException e) {
            // Exception when asked tags doesn't exist
        }
    }

    public void liveStarted() {}
    public void liveStopped() {}
}
```